



IMPORTANT PRODUCT INFORMATION FOR BIOMET 3i BELLATEK® COMPETITIVE CONNECTION PATIENT SPECIFIC ABUTMENTS

This document applies to BIOMET 3i BellaTek Competitive Connection Patient Specific Abutments. For additional product information, please refer to the accompanying package insert or Restorative Manual. The following compatibility chart is specific to BellaTek Competitive Connection Patient Specific Abutments only.

BellaTek Competitive Connection Patient Specific Abutments*

Manufacturer	Implant	Platform (mm)	Abutment Part Numbers
Nobel Biocare®	Internal Connection NobelReplace® Implant	3.5	NRLDAT _x NRLDAN _x
		4.3	
		5	
	Internal Connection NobelActive® Implant	6	NALDAT _x NALDAN _x
		3.5	
		4.3	
Straumann®	Straumann® Bone Level Implant	5.0	SBLDAT _x SBLDAN _x
		3.3	
		4.1	
		4.8	

*NOTE: The maximum angle of the BellaTek Competitive Connection Patient Specific Abutments should not exceed 30°.



BIOMET 3i, LLC
 4555 Riverside Drive
 Palm Beach Gardens, FL 33410
 1-800-342-5454
 Outside The U.S.: +1-561-776-6700
 Fax: +1-561-776-1272
www.biomet3i.com

BIOMET 3i Dental Iberica, S.L.
 WTC Almeda Park, Ed. 1, Planta 1a
 Pl. de la Pau s/n
 08940 – Cornellà de Llobregat
 (Barcelona) Spain
 Phone: +34 934 705 500
 Fax: +34 933 717 849



: Non- sterile	: Lot number
: Catalogue number	: Consult Instructions for use www.ifu.biomet3i.com
: Do not reuse	: Caution: Law Prohibits dispensing without prescription
: European Authorized Representative	: Date of Manufacture



P-BELSUP

P-BELSUP Rev. A
 01/2014